Operating Instructions

863 Valley View Road, Eighty Four, PA 15330
724-941-9701 • skcinc.com

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Checking BioStage/Kit Contents

Use the table below to verify that you received all items associated with the Cat. No. ordered. If you are missing items, contact SKC at 800-752-8472 (U.S. only) or 724-941-9701.

If You Ordered Cat. No.	Your Package Should Contain
225-9611	BioStage [‡] single-stage viable cascade impactor
228-9530K	BioStage Pump Kit includes BioStage [‡] Sampler, QuickTake 30 ^Δ pump with Li-Ion battery, AC charger/adaptor (100-240 V), mounting bracket with inlet adapter, calibration adapter, cassette/tubing adapters, rotameter, and deluxe carry case

[‡] Requires microbiological media supplied by analytical laboratories. For lab list, go to www.skcinc.com/lablocations.

^Δ QuickTake 30 pump is **not** CE marked.

Required Equipment/Media

- Appropriate agar in a 15 x 90 to 100-mm agar plate**
Plastic or glass agar plates can be used. Verify with your agar supplier that the plates contain the proper volume of agar to achieve the appropriate agar height to maintain impactor cut-points. Store agar plates as directed by the supplier.

Suggested Media

For Bacteria: Tryptic Soy Agar (TSA) or Blood Agar Plates (BAP)

For Mold: Potato Dextrose Agar (PDA), Malt Extract Agar (MEA), Dichloran Glycerol 18 Agar (DG-18), or Corn Meal Agar (CMA)

- High flow pump, such as QuickTake 30 Sample Pump Cat. No. 228-9530**
- 1/4-inch ID tubing**
- 3/8-inch ID tubing**
- Luer adapters for 1/4 and 3/8-inch tubing**

For information on laboratories that can provide agar plates and analyze samples, search on Lab at www.skcinc.com.

Clean the BioStage before the first use and between subsequent uses. See Clean and Sterilize.

Sanitize hands or wear disposable gloves. Sanitize hands and impactor any time contamination from handling is possible. Do not touch holes in the jet classification stage.

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
4. Gently screw the jet classification stage back onto the base plate. Align and press the inlet cone onto the jet classification stage until a secure seal is established.

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Pump inlet adapter

QuickTake 30 pump and BioStage in calibration train

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Accessories/Replacement Parts

Description	Cat. No.
QuickTake 30 Sample Pump, Rotameter, and Charger	228-9530A
100-240 V	
Calibration Adapter for BioStage, allows tubing to connect to BioStage inlet	P33100
Mounting Bracket for QuickTake 30, holds BioStage in place on pump during sampling, includes thumbscrews and inlet adapter	228-9531
Tripod Stand for BioStage Impactor	228-506
Replacement O-ring for Base Plate, ea	P32287
Replacement O-ring for Inlet Cone, pk/2	P31893

SKC Limited Warranty and Return Policy

SKC products are subject to the SKC Limited Warranty and Return Policy, which provides SKC's sole liability and the buyer's exclusive remedy. To view the complete SKC Limited Warranty and Return Policy, go to skcinc.com/warranty.

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Macher, J., (ed.) *Bioaerosols: Assessment and Control*, ACGIH, 1999

Samimi, B. and Shufutinsky, A., “Comparison of the Thermo-Andersen N6, the Aerotech A6, the SKC BioStage, and the SKC Micromedia Viable Samplers in Collecting Airborne Fungal Spores,” *AHce 2005*, Anaheim, CA, Final Program, p. 43

Yao, M. and Mainelis, G., “Analysis of Portable Impactor Performance for Enumeration of Viable Bioaerosols,” *Journal of Occupational and Environmental Hygiene*, 4:7, July 2007, pp. 514-524, DOI: 10.1080/15459620701407388

Dobranic, J., “Superbugs in Our Communities - An Introduction for the IEQ Professional,” *EMSL, IAQA Conference*, 2008, Tampa, FL

skcinc.com

Form 40083 Rev 220712